
# ***Rapid New Coronavirus (COVID19) IgG/IgM Combo Test***

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**IMMUNOCHROMATOGRAPHIC ASSAY FOR THE QUALITATIVE ASSESSMENT OF IgM and IgM ANTI-nCORONAVIRUS 19 (COVID19) ANTIBODIES IN WHOLE BLOOD, SERUM AND PLASMA SAMPLES.**

**INTENDED USE**

**Rapid Coronavirus 19 (COVID19) IgG/IgM combo Test** is an *in vitro* qualitative immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to novel *Coronavirus 19 (COVID19) virus* in human whole blood, serum and plasma specimen.

**INTRODUCTION**

Coronavirus (CoV) belongs to the genus Nestovirus, Coronaviridae, and is divided into three genera:  $\alpha$ ,  $\beta$ , and  $\gamma$ . The genus  $\alpha$  and  $\beta$  are only pathogenic to mammals. The genus  $\gamma$  mainly causes bird infections. CoV is mainly transmitted through direct contact with secretions or through aerosols and droplets. There is also evidence that it can be transmitted through the fecal-oral route. So far, there are 7 types of human coronavirus (HCoV) that cause human respiratory diseases: HCoV-229E, HCoV-OC43, SARS-CoV, HCoV-NL63, HCoV-HKU1, MERS-CoV and new coronaviruses (2019) , Is an important pathogen of human respiratory infections. Among them, the new coronavirus (2019) was discovered due to Wuhan virus pneumonia cases in 2019. The clinical manifestations are systemic symptoms such as fever and fatigue, accompanied by dry cough and dyspnea, etc., which can rapidly develop into severe pneumonia, respiratory failure, and acute breathing. Distress syndrome, septic shock, multiple organ failure, severe acid-base metabolism disorders, etc. are even life-threatening.

**TEST PRINCIPLE**

**Rapid Coronavirus 19 (COVID19) IgG/IgM combo Test** utilizes the principle of Immunochromatography. Mouse anti-human IgM and human IgG antibodies are immobilized on the nitrocellulose membrane respectively, as two individual test lines (IgM line and IgG line) in the test window of the test device. As the test sample flows through the membrane within the test device, the colored-Coronavirus 19 (COVID19) specific recombinant antigen-colloidal gold conjugate complexes with specific antibodies (IgM and/or IgG) of Coronavirus 19 (COVID19) virus, if present in the sample. This complex moves further on the membrane to the test region where it is captured by the anti-human IgM and/or human IgG antibodies coated on the membrane leading to formation of a colored band, which indicates a positive test results. Absence of this colored band in the test window indicates a negative test result. A built-in control line will appear in the test window when the test has been performed properly, regardless if positive or negative.

**REAGENTS AND MATERIALS PROVIDED**

Each kit contains:

- |   |            |
|---|------------|
| 1. <b>Rapid Coronavirus 19 (COVID19) IgG/IgM Combo Test</b> | 25 devices |
| 2. Sample buffer  | 1 vial     |
| 3. Instruction for use                                      | 1 piece    |

**MATERIALS NOT PROVIDED**

1. Specimen collection container.
2. Pipetter
3. Timer

**STORAGE AND STABILITY**

The sealed pouches in the test kit may be stored between 4-30°C for the duration of the shelf life as indicated on the pouch.

**PRECAUTIONS**

1. This kit is for *in vitro* diagnostic use only.
2. This kit is for professional use only.
3. Read the instructions for use carefully before performing the test.

**RTVI039 Rapid New Coronavirus 19 (COVID19) IgG/IgM combo test**

Info: sales@immunospark.com

4. This product does not contain any human source materials. Do not use kit after the expiration date.
5. Handle all specimens as potentially infectious.
6. Follow standard Lab procedure and biosafety guidelines. Dispose specimens after autoclaving the at 121°C for at least 20 minutes. Alternatively, treat with 0,5% sodium Hypochlorite for 1-2 hours.
7. Do not pipette reagent by mouth and no smoking or eating while performing assays.
8. Wear gloves during the whole procedure.

**SPECIMEN COLLECTION AND PREPARATION**

1. No prior special preparation of the patient is required before sample collection.
2. If testing cannot be performed immediately, serum/plasma may be stored at 2-8°C up to 3 days. For long-term storage at -20°C for 3 months or -70°C for longer period. Samples with a anticoagulant like EDTA or Heparin or Oxalate may stored at 2-8°C up to 3 days. Blood samples should not be frozen.
3. Repeated freezing and thawing of the specimen should be avoided
4. Heat inactivation of specimens should be avoided.
5. Avoid to use haemolysed, clotted, contaminated, lipamic and viscous/turbid specimen.
6. Specimen containing precipitates must be centrifuged and the clear supernatant only used.
7. Do not inactivate the sample by heating.

**PROCEDURE**

1. Bring the kit components to room temperature before testing.
2. Open the pouch and remove the card. Once opened, the test must be used immediately.
3. Label the test card with the patient identity.
4. **For Serum or Plasma Specimens:**

Using the provided 10uL disposable pipette, draw the specimen up to the Fill Line, and transfer 1 drop of serum/plasma (approximately 5uL),to the specimen well of the test device, then add 1 drop of buffer and start the timer.


**For Whole Blood (Venipuncture/Fingerstick) Specimens:**

Using the provided 10uL disposable pipette, and transfer 2 drops of whole blood (approximately 20µL) to the specimen well of the test device, then add 1 drop of buffer and start the timer.

5. Wait for the colored line(s) to appear. Read results at **15 minutes. Results interpreted after 20 minute may not be accurate.**

**INTERPRETATION OF RESULTS**

- 1) **IgG Positive:** The presence of two color bands, “C” and “G”, indicates a positive result for IgG.
- 2) **IgM Positive:** The presence of two color bands, “C” and “M” indicates a positive result for IgM.
- 3) **IgG and IgM Positive:** The presence of three color bands indicates a positive result for IgG and IgM.
- 4) **Negative:** The presence of only one band, “C” ,within the result window indicates a negative result
- 5) **Invalid:** Test is invalid if the control line does not appear. If this occurs, the test should be repeated using a new device. If the problem persists refer to your local distributor.

**QUALITY CONTROL**

1. The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials which is not provided with this test kit may be commercially available.

**RTVI039 Rapid New Coronavirus 19 (COVID19) IgG/IgM combo test**

Info: sales@immunospark.com

**LIMITATIONS**

1. The test is qualitative detection of anti-CoV 19 (COVID19) antibody in serum, plasma or whole blood.
2. A definitive clinical diagnosis should not be based on the result of a single test.

**PERFORMANCE CHARACTERISTICS**

**Diagnostic Sensitivity and Specificity**

**Rapid New Coronavirus 19 (COVID19) IgG/IgM combo Test** has been evaluated using a total amount of 400 samples from patients with suspect of infection (symptomatics) verified with a Real-PCR assay:

		Realtime-PCR				
		Positive	Negative	Total		
<b>RTVI039 nCOV19 IgM Test</b>					<b>Diagnostic sensitivity</b>	93,46%
	<b>Positive</b>	100	2	102	<b>Diagnostic specificity</b>	99,32%
	<b>Negative</b>	7	291	298	<b>Overall agreement</b>	97,75%
	<b>Total</b>	107	293	<b>400</b>	<b>Confidence interval</b>	95,00%

		Realtime-PCR				
		Positive	Negative	Total		
<b>RTVI039 nCOV19 IgG Test</b>					<b>Diagnostic sensitivity</b>	98,15%
	<b>Positive</b>	106	4	110	<b>Diagnostic specificity</b>	98,63%
	<b>Negative</b>	2	288	290	<b>Overall agreement</b>	98,50%
	<b>Total</b>	108	292	<b>400</b>	<b>Confidence interval</b>	95,00%

**Diagnostic accuracy**

**Healthy patients (verified as negative with the swab):**

**Rapid New Coronavirus 19 (COVID19) IgG/IgM combo Test** has been compared with a RT-PCR assay on 100 samples from healthy patients:

Results	RTVI039 IgG	RTVI039 IgM	RT-PCR
<b>Positive</b>	0	0	0
<b>Negative</b>	100	100	100%
<b>Accuracy</b>	100%	100%	100%

**Infected patients (verified as positive with the swab):**

**Rapid New Coronavirus 19 (COVID19) IgG/IgM combo Test** has been compared with a RT-PCR assay on 200 samples from patients with confirmed infection:

Results	RTVI039 IgG	RTVI039 IgM	RT-PCR
<b>Positive</b>	190	183	190
<b>Negative</b>	10	17	10
<b>Accuracy</b>	95%	91,5%	95%

**Interferences**

Negative samples have been spilled with the following substances. No interferences have been found at the following concentrations:

<b>Hemoglobine</b>	5 g/L	<b>Tryglicerides</b>	3.89 mmol/L
<b>Bilirubine</b>	20 mg/dL	<b>Cholesterol</b>	6.10 mmol/L

**Repeatability**

A pool of positive samples to IgG and to IgM has been submitted to 2 sequential dilutions (1/50 and 1/100), additionally a pool of confirmed negative samples has been considered.

The 3 pools have been tested in 10 replicates with **RTVI039**.

**Reproducibility**

The same 3 pools of samples have been tested in 10 replicates using 3 different lots of **RTVI039**.

The results reproduced 99,99%.

**BIBLIOGRAFIA/REFERENCES**

1. World Health Organization (WHO), WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China, Beijing: WHO; 9 Jan 2020.
2. World Health Organization (WHO), Report of the WHO-China Joint Mission on Coronavirus Disease 2019 (COVID-19), WHO; 28 Feb 2020.
3. Fan Y, Zhao K, Shi ZL, Zhou P (March 2019). "Bat Coronaviruses in China". *Viruses*. 11 (3): 210. doi:10.3390/v11030210. PMC 6466186. PMID 30832341.
4. De Groot RJ, Baker SC, Baric R, Enjuanes L, Gorbalenya AE, Holmes KV, Perlman S, Poon L, Rottier PJ, Talbot PJ, Woo PC, Ziebuhr J (2011). "Family *Coronaviridae*". In King AM, Lefkowitz E, Adams MJ, Carstens EB, International Committee on Taxonomy of Viruses, International Union of Microbiological Societies. Virology Division (eds.). *Ninth Report of the International Committee on Taxonomy of Viruses*. Oxford: Elsevier. pp. 806–28. ISBN 978-0-12-384684-6.
5. Kahn, Jeffrey; McIntosh, Kenneth (November 2005), "History and recent advances in coronavirus discovery", *Pediatric Infectious Disease Journal*, 24 (11): s223–s227, doi:10.1097/01.inf.0000188166.17324.60, archived from the original on 2020-02-05.
6. World Health Organization (28 January 2020). "Clinical management of severe acute respiratory infection when novel coronavirus (2019-nCoV) infection is suspected"

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HEADQUARTER: Via lucrino 35 – 00199 – Rome – Italy

 e-mail: [sales@immunospark.com](mailto:sales@immunospark.com)

 R&D and MANUFACTURING: Via del mare 187, 00071, Pomezia – Rome – Italy;  
 phone + 39 06 64012332 - Fax + 39 06 21113989; e-mail: [prod.pomezia@immunospark.com](mailto:prod.pomezia@immunospark.com)
[www.immunospark.com](http://www.immunospark.com)